Benefits of early in vitro screening for seizure liability in problem solving and decision making

K L Rockley¹; R A Roberts¹; M J Morton¹

¹ApconiX, Alderley Park, Alderley Edge, Cheshire, UK



Seizure liability remains a significant cause of attrition throughout drug development both in pre-clinical and clinical studies. This emphasizes the need for improved methodologies to detect seizure liability prior to in vivo toxicology studies, ideally with reduced reliance on animals and better translation to humans. Much like the Comprehensive in vitro Proarrhythmia Assay (CiPA) which is now widely accepted for early assessment of cardiovascular safety, we have developed an approach utilizing hiPSC-neuronal cell microelectrode array (MEA) and ion channel screening for early seizure prediction. These studies provide mechanistic information, and are useful for early de-risking, and support optimal drug design.

EXAMPLES DEMONSTRATING THE UTILITY OF IN VITRO SEIZURE LIABILITY ASSAYS

Known Seizurogenic (1)compounds identified correctly with high predictivity

(2) **De-risk and prioritize a** chemical series at single concentration

(3) Good correlations between *in vitro - in vivo* responses with novel compounds

(4) Assess human relevance of convulsions observed in nonclinical studies

(1) PREDICITIVITY OF MEA ASSAY

- ✓ 14/16 known seizurogenic compounds were identified correctly (87%)
- Correlations between in vitro responses and clinical exposures
- ✓ E.g. amoxapine An antidepressant where 8.75% ADRs are convulsions

				Amoxapine				Baseline	Amoxapine
	MEA parameters	0.1µM	0.3µM	1μM	3μΜ	10μΜ			
Cnike	Mean firing rate (Hz)	NC	NC	↑	$\uparrow\uparrow$	$\downarrow \downarrow \downarrow \downarrow$			
зріке	Interspike Interval (ISI) Coefficient of variation - Avg	NC	NC	NC	\rightarrow	$\downarrow \downarrow \downarrow \downarrow$	10uM	and the same size and same size after all same	and a set of second
	Burst frequency - Avg (Hz)	NC	NC	ተተ	$\uparrow \uparrow \uparrow$	$\downarrow \downarrow \downarrow \downarrow$	TOPIN		
Durata	Burst duration - Avg (sec)	NC	NC	\downarrow	$\rightarrow \rightarrow$	$\downarrow \downarrow \downarrow \downarrow$			
Bursts	Number of spikes per burst - Avg	NC	NC	\downarrow	$\downarrow\downarrow$	$\downarrow \downarrow \downarrow \downarrow$		8 50 130 191 Time (sec)	0 k0 100 Time (sec)
	Mean ISI within burst - Avg (sec)	NC	NC	NC	$\downarrow\downarrow$	UNC			
	Network burst frequency	NC	NC	ተተ	$\uparrow\uparrow\uparrow$	N/A			AUDUCAL MARKING RECORD OF THE
	Network burst duration - Avg (sec)	NC	NC	\downarrow	\rightarrow	N/A	ЗμМ		
Network Bursts	Number of spikes per network burst - Avg	NC	NC	\downarrow	$\downarrow \downarrow \downarrow \downarrow$	N/A			
	Mean ISI within network burst - Avg (sec)	NC	NC	NC	$\uparrow\uparrow$	N/A		111111111	
	Network IBI Coefficient of variation - Avg	NC	NC	NC	$\downarrow\downarrow$	UNC		9 00 100 100	Time (sec)
Synchrony	Area under normalized cross-correlation	NC	NC	NC	NC	$\downarrow \downarrow \downarrow \downarrow$			
NC no cl (±109 one arrow 0 – 2 two arrows 30 –	Therapeutic and toxic doses taken from Schulz et al., 2 %) and toxic blood concentrations of more than 1100 dru 9% change	2020 – Revisited: Th gs and other xenob	erapeutic Th iotics	nerapeutic ran 0.57 - 1.91μΝ	ge	Toxic dose 9.56μM	1μM		



(2) DE-RISKING A CHEMICAL SERIES

Testing a chemical series revealed one with low seizure risk (compound A)

This compound had distinct structural features.



(3) IN VITRO – IN VIVO TRANSLATION

	NBI compound B						
	MEA parameters	1μΜ	3µM	10µM	30µM	100μΜ	
Creike	Weighted mean firing rate (Hz)	NC	NC	\downarrow	\downarrow	$\downarrow\downarrow$	 Minimal changes at
Бріке	Interspike Interval (ISI) Coefficient of variation - Avg	NC	NC	NC	NC	NC	
	Burst frequency - Avg (Hz) *	\uparrow	\uparrow	\uparrow	\uparrow	\uparrow	lower doses
Bursts	Burst duration - Avg (sec)	NC	\downarrow	\downarrow	\downarrow	$\checkmark \checkmark$	
	Number of spikes per burst - Avg	NC	J.	\downarrow	\downarrow	*	

$(\mathbf{4})$ HUMAN RELEVANCE OF CONVULSIONS IN DIFFERENT SPECIES

Nonclinical testing of a compound caused convulsions only in dogs On the MEA assay only the dog-specific metabolite caused a seizurogenic phenotype



and 2.54 μ M (total) (1000mg/kg, PO)

 \checkmark Only the dog metabolite hit an ion channel on the panel providing mechanistic insight into the possible cause of seizure

		Parent drug							
	MEA parameters	1μΜ	3μΜ	10µM	30µM	100μΜ	•		
Spike	Weighted mean firing rate (Hz)	NC	NC	\uparrow	NC	NC			
	Interspike Interval (ISI) Coefficient of variation - Avg	\uparrow	NC	NC	NC	\uparrow			
Dourse	Burst frequency - Avg (Hz)	NC	NC	NC	NC	NC			
	Burst duration - Avg (sec)	NC	NC	\uparrow	\uparrow	\uparrow			
DUISIS	Number of spikes per burst - Avg	NC	NC	个个	\uparrow	\uparrow			
	Mean ISI within burst - Avg (sec)	NC	NC	\downarrow	\checkmark	NC	•		
Network Bursts	Network burst frequency	NC	NC	NC	NC	\rightarrow			
	Network burst duration - Avg (sec)	NC	NC	NC	NC	NC			
	Number of spikes per network burst - Avg	NC	NC	NC	\uparrow	NC			
	Mean ISI within network burst - Avg (sec)	NC	NC	NC	NC	NC			
	Network IBI Coefficient of variation - Avg	\rightarrow	\uparrow	\uparrow	\downarrow	UNC			
Synchrony	Area under normalized cross-correlation	NC	NC	NC	NC	NC			

Inactive up to 100µM

No hits on ion channel panel

Human and rodent metabolite

MEA parameters	1μΜ	3μΜ	10μΜ	30µM	100µM
Weighted mean firing rate (Hz)	NC	NC	NC	NC	\rightarrow
Interspike Interval (ISI) Coefficient of variation - Avg	NC	NC	NC	Λ 30μΜ NC NC NC √ NC √ NC NC	NC
Burst frequency - Avg (Hz)	NC	NC	\rightarrow	\rightarrow	\rightarrow
Burst duration - Avg (sec)	NC	NC	NC	NC	NC
Number of spikes per burst - Avg	NC	NC	NC	NC	NC
Mean ISI within burst - Avg (sec)	NC	NC	NC	M 30μM 30 C NC 0 C NC 0 Λ ↓ 0 Λ ↓ 0 Λ ↓ 0 Λ ↓ 0 Λ ↓ 0 Λ ↓ 0 Λ NC 0 Λ ΝC 0 Λ Λ 0 Λ ΝC 0 Λ Λ 0 Λ Λ 0 Λ Λ 0 Λ Λ 0 Λ Λ 0 Λ Λ 0 Λ Λ 0 Λ Λ 0 Λ Λ 0 </td <td>\uparrow</td>	\uparrow
Network burst frequency	NC	NC	\downarrow	NC	\rightarrow
Network burst duration - Avg (sec)	NC	NC	NC	NC	NC
Number of spikes per network burst - Avg	NC	NC	NC	NC	NC
Mean ISI within network burst - Avg (sec)	NC	NC	NC	\uparrow	\uparrow
Network IBI Coefficient of variation - Avg	UNC	UNC	UNC	UNC	UNC
Area under normalized cross-correlation	NC	NC	NC	NC	\downarrow
	MEA parametersWeighted mean firing rate (Hz)Interspike Interval (ISI) Coefficient of variation - AvgBurst frequency - Avg (Hz)Burst duration - Avg (sec)Number of spikes per burst - AvgMean ISI within burst - Avg (sec)Network burst frequencyNetwork burst duration - Avg (sec)Number of spikes per network burst - AvgMean ISI within network burst - Avg (sec)Number of spikes per network burst - AvgMean ISI within network burst - Avg (sec)Network IBI Coefficient of variation - AvgArea under normalized cross-correlation	MEA parameters1µMWeighted mean firing rate (Hz)NCInterspike Interval (ISI) Coefficient of variation - AvgNCBurst frequency - Avg (Hz)NCBurst duration - Avg (sec)NCNumber of spikes per burst - AvgNCMean ISI within burst - Avg (sec)NCNetwork burst frequencyNCNetwork burst frequencyNCNumber of spikes per network burst - Avg (sec)NCNetwork burst duration - Avg (sec)NCNumber of spikes per network burst - AvgNCNetwork lurst duration - Avg (sec)NCNetwork IBI Coefficient of variation - AvgUNCArea under normalized cross-correlationNC	MEA parameters1μM3μMWeighted mean firing rate (Hz)NCNCInterspike Interval (ISI) Coefficient of variation - AvgNCNCBurst frequency - Avg (Hz)NCNCBurst duration - Avg (sec)NCNCNumber of spikes per burst - AvgNCNCMean ISI within burst - Avg (sec)NCNCNetwork burst frequencyNCNCNetwork burst duration - Avg (sec)NCNCNetwork burst frequencyNCNCNetwork burst duration - Avg (sec)NCNCNumber of spikes per network burst - AvgNCNCNumber of spikes per network burst - Avg (sec)NCNCNumber of spikes per network burst - Avg (sec)NCNCNetwork IBI Coefficient of variation - AvgUNCUNCArea under normalized cross-correlationNCNC	MEA parameters1μM3μM10μMWeighted mean firing rate (Hz)NCNCNCInterspike Interval (ISI) Coefficient of variation - AvgNCNCNCBurst frequency - Avg (Hz)NCNC√Burst duration - Avg (sec)NCNCNCNumber of spikes per burst - AvgNCNCNCMean ISI within burst - Avg (sec)NCNCNCNetwork burst frequencyNCNCNCNetwork burst frequencyNCNCNCNumber of spikes per network burst - AvgNCNCNCNetwork burst duration - Avg (sec)NCNCNCNumber of spikes per network burst - AvgNCNCNCNumber of spikes per network burst - AvgNCNCNCNetwork IBI Coefficient of variation - AvgUNCUNCUNCArea under normalized cross-correlationNCNCNC	MEA parameters1μM3μM10μM30μMWeighted mean firing rate (Hz)NCNCNCNCInterspike Interval (ISI) Coefficient of variation - AvgNCNCNCNCBurst frequency - Avg (Hz)NCNCNCVVBurst duration - Avg (sec)NCNCNCNCNCNumber of spikes per burst - AvgNCNCNCNCNCMean ISI within burst - Avg (sec)NCNCNCNCNCNetwork burst duration - Avg (sec)NCNCNCNCNCNumber of spikes per network burst - AvgNCNCNCNCNumber of spikes per network burst - Avg (sec)NCNCNCNCNetwork IBI Coefficient of variation - AvgUNCUNCUNCUNCArea under normalized cross-correlationNCNCNCNCNC

Inactive up to 100µM

No hits on ion channel panel

Dog metabolite

MEA parameters		1μΜ	3μΜ	10µM	30µM	100µM
Calles	Weighted mean firing rate (Hz)	NC	NC	NC	\rightarrow	NC*
зріке	Interspike Interval (ISI) Coefficient of variation - Avg	NC	NC	NC	NC	\rightarrow
	Burst frequency - Avg (Hz)	NC	NC	\uparrow	\uparrow	$\uparrow\uparrow$

Clear seizurogenic phenotype at high



DISCUSSION AND CONCLUSIONS

- Currently, the majority of seizure detection for de-risking is done in animal models, with limited success and questionable translation to humans.
- Our approach has been validated using known seizurogenic compounds and shows good translation to clinical exposures (Rockley et al., 2023).
- These assays are useful for early de-risking and problem solving later in drug development.
- Collectively, these studies demonstrate the utility of this approach for early seizure prediction to provide mechanistic information, early de-risking, and support optimal drug design using human in vitro models.